

***Crystallization Demonstration***  
***Phenyl Salicylate Crystals***  
**Teacher Resource**

**Grade Level**      **6-12**

**Objectives**

1. The student will develop an understanding of temperature and the melting points of chemicals.
2. The student will employ simple equipment and tools to gather data.
3. The student will demonstrate the use of 21<sup>st</sup> century technology with a digital microscope to locate, evaluate, and collect information.

**National Standards**

<a href="#">NS.K-4.1</a> ; <a href="#">NS.5-8.1</a> ; <a href="#">NS.9-12.1</a>	Science as Inquiry
<a href="#">NS.K-4.2</a> ; <a href="#">NS.5-8.2</a> ; <a href="#">NS.9-12.2</a>	Physical Science
<a href="#">NS.K-4.5</a> ; <a href="#">NS.5-8.5</a> ; <a href="#">NS.9-12.5</a>	Science & Technology
<a href="#">NT.K-12.1</a>	Basic Operations and Concepts
<a href="#">NT.K-12.3</a>	Technology Productivity Tools
<a href="#">NT.K-12.6</a>	Technology Problem-Solving & Decision-Making Tools

**Teacher Background Information**

This demonstration is perfect to use as a demonstration to a unit on crystallization. Crystallization occurs through a natural process where solid crystals precipitate from either a solution, by melting, or by a gas. This demonstration will show crystallization occurring by heating phenyl salicylate in its solid form around 43°C, so the process is through cooling. When it begins to lower in temperature from the melting point, crystals will begin to form.

**MATERIALS**

Digital Microscope	Glass Microscope Slides
Personal computer	Hot Plate
Interactive White Board/Projector	Forceps
Phenyl Salicylate	Safety goggles

## PROCEDURE

1. Wear safety goggles at all times.
2. Sprinkle a small amount of phenyl salicylate on a glass slide.
3. Heat on a hot plate until the phenyl salicylate powder starts to melt.
4. Immediately remove the slide from the hot plate with the forceps.

*Note: place 2 washers on the hot plate with the slide on top of the washers...this will be easier to remove the slide from the hot plate.*

5. Place the slide underneath the digital microscope and focus.
6. Start recording the recrystallization of phenyl salicylate under the microscope.

## DISCUSSION

Connect this demonstration to geologic processes with your class about crystallization, the differences between minerals cooling rates, and the orderly framework of how minerals form (these ideas can be seen in the demonstration).