


Discovery Scope Kit

ASSEMBLY & USE

EXPLORER DISCOVERY SCOPE KIT

For use with Video Flex® 2000 series



ASSEMBLY INSTRUCTIONS:


EXPLORER DISCOVERY SCOPE KIT (VFDS2000)

1. Attach Explorer clamp onto the neck of the Video Flex®, spaced approximately 1" down from camera head. Gently tighten the wing-nut down on the tweezer clamp.
2. Attach the Telescoping Tube to the other end of explorer clamp as shown above. Gently tighten the wing-nut on the 2nd tweezer-clamp.
3. Select which attachment you would like to use (either the tweezer-clamp for 35mm slides and other grippable items, or the specimen "Bug" box holder), then slide the attachment shaft through the opening located at the end of the telescoping tube
4. Adjust the telescoping tube and attachment holders to position the object directly in front of the Video Flex® camera lens. You can quickly and easily interchange the attachment holders to suit your needs.

You are only limited by your imagination!

DISCOVERY SCOPE KIT

For use with all Video Flex® except 2000 series



DISCOVERY SCOPE KIT (VFDS)

1. Attach the Discovery Scope Clamp onto the neck of the Video Flex®, space approximately 1/2" down from camera head. Simply push clamp onto neck.
2. Attach the Telescoping Tube to the other clamp opening as shown above.
3. Select which attachment you would like to use either the tweezer-clamp for the 35mm slides and other grippable items, or the specimen "Bug" box holder, then slide the attachment shaft through the opening located at the end of the telescoping tube.
4. Adjust the telescoping tube and attachment holders to position the object directly in front of the Video Flex® camera lens. Voila-that's it! You can quickly and easily interchange the attachment holders to suit your needs.

You are only limited by your imagination!



5615 RAYTOWN ROAD-KANSAS CITY, MO 64133

816.535.4787-FAX 816.358.5072

E-mail: info@ken-a-vision.com <http://www.ken-a-vision.com>

Learn more about the complete line of Video Flex® multi-purpose flexible video microscopes, Microprojectors, and Microscopes-call 1-800-334-8407 or explore our website at: www.ken-a-vision.com